
E11.OGL.002

Calomel Electrode

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E11.OGL.002

Calomel Electrode

Introduction

The E11.OGL.002 is a Calomel Reference Electrode fitted with a 14,5/23 sleeve.

Preparation for Measurement

1. Remove the protection cap from the electrode and any seals covering the filling hole and the porous plug (pin).
2. Before starting a measurement, remove the clip which closes the electrode filling hole. Remember to replace the clip at the end of measurements.
3. Check the level of the filling solution. It should be approximately 0.5 cm below the filling hole.
4. In order to remove air bubbles trapped inside the electrode, shake the electrode holding it at its head with the porous pin down.

Maintenance

1. Electrode contamination is a major cause of faulty measurements.
2. The electrode should be rinsed with distilled water after measurements.
3. Check frequently the level of filling solution.
4. In case of deposits which cover the electrode, clean the electrode with:
 - a solution of acid (0.1M HCl, 0.1M HNO₃): mineral salt deposits, etc...
 - Normal Cleaning Solution : greasy or oily deposits.

The porous plug of the electrode can be cleaned using a fine abrasive paper.

Storage

Between measurements: leave it in a Saturated KCl Solution.

Several days: seal the filling hole with paraffin film or with the electrode clip and immerse the electrode in a Saturated KCl Solution using the accompanying cap.

Accessories

- ✓ CAD015RV1
Electrode Cable, with BNC plug, 1m
- ✓ CAD011RV1
Electrode Cable, with BNC plug, 0,16m
- ✓ CAD032RV1
Electrode Cable, with Banana \varnothing 4mm plug, 1m
- ✓ Normal Cleaning solution 150 ml
- ✓ Normal Storage solution 150 ml
- ✓ Saturated KCL Solution 150 ml

Specifications

Temperature range: -10 to 50°C

N.B.: In case of contact with mercury, wash exposed areas immediately with plenty of water.

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Origalys ElectroChem SAS