
## E10.OGL.026

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# Reference Electrode

### Introduction

The E10.OGL.026 is a Hg/HgO Reference Electrode, ideal for use in alkaline solutions.

### Preparation for Measurement

1. Remove the protection cap from the electrode and any seals covering the filling hole and the porous plug (pin).
2. Before starting a measurement, remove the clip which closes the electrode filling hole. Remember to replace the clip at the end of measurements.
3. Check the level of the filling solution. It should be approximately 0.5 cm below the filling hole.
4. For tests and on delivery, the electrodes are filled with 0.1 M potassium hydroxide. However, we highly recommend that you replace this solution with a solution which has a better stability, e.g. KOH 1M. If carbonates are formed in the KOH solution, the reference potential will drift. Other concentrations or other hydroxides may therefore be used.

Note: extreme care must be taken when working with strong concentrations of potassium hydroxide as they are very corrosive. In case of contact with the skin, wash immediately with plenty of water.

5. In order to remove air bubbles trapped inside the electrode, shake the electrode holding it at its head with the porous pin down.

### Maintenance

1. Electrode contamination is a major cause of faulty measurements.
2. The electrode should be rinsed with distilled water after measurements.
3. Check frequently the level of filling solution.

4. In case of deposits which cover the electrode, clean the electrode with:

- a solution of acid (0.1M HCl, 0.1M HNO<sub>3</sub>): mineral salt deposits, etc...
- Normal Cleaning Solution : greasy or oily deposits.

### Storage

**Between measurements:** leave it in a potassium hydroxide solution.

**Several days:** cover the porous pin or fibres with paraffin film or with the electrode clip and immerse the electrode in potassium hydroxide solution.

**Long-term storage:** it is recommended to empty the electrode and refill it with a fresh solution of KOH 0.1M. Do not forget to fill the protection cap with paraffin film and place the electrode in its box and store at room temperature.

### Accessories

- ✓ CAD015RV1  
Electrode Cable, with BNC plug, 1m
- ✓ CAD011RV1  
Electrode Cable, with BNC plug, 0,16m
- ✓ CAD032RV1  
Electrode Cable, with Banana Ø 4mm plug, 1m
- ✓ Normal Cleaning solution 150 ml
- ✓ Normal Storage solution 150 ml
- ✓ Saturated KCL Solution 150 ml

### Specifications

Temperature range: -5 to 50°C

Note: the reference potential at 25°C depends on the type of filling solution used:

- For an electrode filled with 1M KOH it is approximately -129 mV,
- For an electrode filled with 0.1M KOH it is approximately -70 mV, with respect to a saturated calomel electrode

**N.B.: In case of contact with mercury, wash exposed areas immediately with plenty of water.**

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