
E11.OGL.014

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Reference Electrode

Introduction

The E11.OGL.014 is an Ag/AgCl Reference Electrode, fitted with a screw cap.

Preparation for Measurement

1. Remove the protection cap from the electrode and any seals covering the filling hole and the porous plug (pin).
2. Before starting a measurement, remove the clip which closes the electrode filling hole. Remember to replace the clip at the end of measurements.
3. Check the level of the filling solution. It should be approximately 0.5 cm below the filling hole. If necessary, refill with 3 M KCl Solution, saturated with AgCl.

If desired, the concentration of KCl can be changed. However, it is advisable to use a high concentration. Remember always to saturate the solution with AgCl. For nonaqueous applications such as measurements in acetic acid. Empty the electrode and fill it up again with a saturated KCl solution in acetic acid. Saturate the solution with AgCl.

4. In order to remove air bubbles trapped inside the electrode, shake the electrode holding it at its head with the porous pin down.

Maintenance

1. Electrode contamination is a major cause of faulty measurements.
2. The electrode should be rinsed with distilled water after measurements.
3. Check frequently the level of filling solution.
4. In case of deposits which cover the electrode, clean the electrode with:
 - a solution of acid (0.1M HCl, 0.1M HNO₃): mineral salt deposits, etc...

- Normal Cleaning Solution : greasy or oily deposits.

The porous plug of the electrode can be cleaned using a fine abrasive paper

Storage

Between measurements: leave it in 3M KCl solution.

Overnight or longer: seal the filling hole with paraffin film or with the electrode clip and fit back in place the protection cap filled with the 3M KCl Solution.

Accessories

- ✓ CAD015RV1
Electrode Cable, with BNC plug, 1m
- ✓ CAD011RV1
Electrode Cable, with BNC plug, 0,16m
- ✓ CAD032RV1
Electrode Cable, with Banana Ø 4mm plug, 1m
- ✓ Normal Cleaning solution 150 ml
- ✓ Normal Storage solution 150 ml
- ✓ 3M KCl saturated with AgCl 150 ml

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Origalys ElectroChem SAS