



APPLICATION NOTE

Determination of diffusion coefficient of electroactive elements through Randles-Sevcik equation by Cyclic Voltammetry

General Electrochemistry AP-GE19

www.originalys.com



Summary:

1 Introduction

2 Parameters

3 Results and discussion

4 References

5 Necessary instruments and materials

Introduction

In electrochemical systems, mass transport of electroactive species through the electrode surface, plays a key role in the shape and intensity of voltametric signals. Diffusion is one of the kinetic factor in mass transport phenomena.

Among the different methods used to investigate diffusion, the Randles-Sevcik equation offers a straightforward way to determine the diffusion coefficient of a redox species from cyclic voltammetry experiments, under conditions of semi-infinite linear diffusion.

When a redox system is electrochemically reversible, the peak current observed in a cyclic voltammogram is directly proportional to the square root of the scan rate.

This relation is described by the Randles-Sevcik equation:

$$i_p = 0,4463nFAC \sqrt{\frac{nFD}{RT}} v_b$$

With:

i_p = Levich current (A)

n = number of electrons

F = Faraday (96,500 C/mol)

A = Surface of Working Electrode (m^2)

C = bulk concentration (mol/ m^3)

D = Ideal gas constant (8,314 J/mol/K)

R = Kinematic viscosity ($10^{-6} m^2/s$)

T = Temperature (K)

v_b = The scan rate (V/s)


The measurements were carried out in a solution containing 0.01 M potassium ferricyanide ($K_3[Fe(CN)_6]$) and 0.5 M potassium chloride (KCl) as supporting electrolyte. The working electrode was a 5 mm diameter platinum disk, offering a well-defined surface area for accurate current measurements and diffusion coefficient determination.

Parameters

This experiment was performed by OrigMaster software. The flow chart and parameters are shown in figure 1. The CV4limits were repeated 20 times increasing the scan rate by 5 mV/s each time (from 5 mV/s to 100 mV/s).

Before sweeping the potential, the OCP of the working electrode was measured for 100 seconds to be sure that the cell assembly is well stable.

Figure 1: Parameters of the Pot. CV 4 limits



TIPS: Enabling the "Rejection 50 Hz" option reduces environmental noise, resulting in cleaner and more readable curves. For further details, please refer to the "Rejection Filter" application note.

Results and Discussion

An overlay of the 20 CV curves made from the previous flowchart is shown in figure 2.




Figure 2: 20 overlaid voltammograms with scan rates varying from 5 mV/s to 100 mV/s

As predicted, it is observed a clear increase in the cathodic peak current with the square root of the scan rate, consistent with diffusion-controlled electron transfer. From the slope of the linear fit, the calculated diffusion coefficient is :

$D = 7.04 \times 10^{-10} \text{ m}^2/\text{s}$, which falls within the expected order of magnitude [1] [2].

These results validate both the experimental setup and the reliability of cyclic voltammetry as a tool for estimating diffusion coefficients. This also reinforces the method's relevance in preliminary studies of electropolishing conditions, where diffusion processes can influence reaction kinetics and surface behaviour.

References

[1] : Klymenko, Oleksiy & Evans, Russell & Hardacre, Christopher & Svir, Irina & Compton, Richard. (2004). Double potential step chronoamperometry at microdisk electrodes: Simulating the case of unequal diffusion coefficients. Journal of electroanalytical chemistry. 571. 211-221. 10.1016/j.jelechem.2004.05.012.

[2] : Nassim Ait Mouheb, Agnès Montillet, Camille Sollicec, Jacques Comiti, Patrick Legentilhomme. Caractérisation des écoulements et des transferts de matière dans les micromélangeurs. CFM 2009 - 19ème Congrès Français de Mécanique, Aug 2009, Marseille, France. ffhal-03391323

Necessary instruments and materials




Code : AR02407

OGS055



Cable and adapter



OrigaMaster 5



S7

Screw Head S7

Code : OGR007

REF: Ag/AgCl Electrode



S7

Screw Head S7

Code : OGV005

AUX: Platinum Electrode



BNC

Code : AR02633CIAL

Platinum Tip ø5 mm
Code: AR02019

WRK: OrigaTrod Lt & tip



Code : AR00484CIAL

Glass cell with support



Code : AR01647CIAL

Polishing kit

- Electrolyte: Ferri Cyanide solution 1×10^{-2} M in KCl 5×10^{-1} M

Have an idea for an Application Note ?

We'd love to hear it!

At OrigalyS, we create tailor-made application notes, adapted to your protocols, methods, and research topics.

Get in touch, let's co-write clear and practical application notes together!



22/07/2025

OrigalyS ElectroChem SAS

555 Chemin du bois
69140 RILLIEUX-LA-PAPE
FRANCE

Tel. +33 (0)9 54 17 56 03

Fax. +33 (0)9 59 17 56 03

contact@origalys.com